

Mobile 2 RPW







Mobile 2 RPW



Chattanooga's Radial Pressure Wave technology is mobile.

The Mobile 2 RPW is compact and light enough to be easily carried while also delivering a host of features and benefits.

When there is need for a powerful yet portable Radial Pressure Wave device in the clinic or on the sports field, the Chattanooga Mobile 2 RPW is the ideal choice. A lightweight and compact design together with a quiet built-in compressor make it perfect for portable use. A user-friendly experience is at the heart of the Mobile 2 RPW's design. The simple, ergonomically laid out controls enable all of the key parameters to be easily changed, while a clear, illuminated LED display removes ambiguity when monitoring pressure, frequency and number of pulses. Pulses can be initiated via a trigger mechanism on the handpiece.



This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- Orthopaedics
- Physiotherapy
- Aesthetic
- Chiropractors

- Rheumatology
- Sport Clubs
- Veterinary

- Physical Medicine
- Acupuncture
- Podiatrists

Indications

- Myofascial pain
- Muscle tendon tissue activation
- Tendon attachment & site illness

- Acupuncture Shockwave Therapy

The Mobile RPW D-Actor Applicator allows quick and efficient uses of applicator transmitters for desired therapeutic effects during treatment.

D-ACTOR® Applicator

Applications:

- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping

· Regulation of muscle tone

Medical massage





Facts and figures

- Weight of the unit: Approx. 9 kg
- Dimensions (L x W x H): 32 x 32 x 21 cm
- Built in compressor
- Single and Continuous shocks
 Frequency D-actor: 1 21 Hz
 Frequency V-actor: 1 31 Hz
- Application pressure: 1 5 Bar
 Highly customisable with precise adjustments of 0.1 increments.
- Oscillating D-Actor
- V-Actor Hand piece option
- USB connectivity



User Interface Features




Hand piece docking port



Ultrasound gel holder



D-ACTOR® Transmitters





Mobile 2 RPW

Ordering Information

| Part Number | Description |
|-------------|--------------|
| 2910 | Mobile 2 RPW |

Standard Accessories (included with Mobile 2 RPW as default

| 2905 | Chattanooga Mobile RPW 2 Unit Only |
|----------|------------------------------------|
| 32000 | R-SW Hand Piece Applicator Set |
| 28178 | Ro40 'Beam' transmitter |
| 29724 | D20-S D-Actor® 20mm Transmitter |
| 13-5182 | Conductor™ Transmission Gel 250 ml |
| 13-00060 | CD User Manual |

Optional Accessories

| 28701 | Revision Kit |
|-------|---|
| 28178 | Ro40 'Beam' transmitter |
| 31800 | 9 Transmitter Set In Case, Accessory |
| 29801 | A6 Point, Transmitter |
| 29802 | T10 Finger, Transmitter |
| 29729 | C15, 15 mm Cerama-X Transmitter, Black |
| 29728 | D115, 15 mm, Deep Impact Transmitter, Black |
| 29726 | F15 Focus-Lens 15 mm Transmitter |
| 29724 | D20-S D-Actor® 20mm Transmitter |
| 28736 | D20-T Transmitter |
| 29725 | D35-S, 35 mm, D-Actor Transmitter, Black |
| 22375 | Handpiece Cushion, Mobile 2 RPW |
| 22376 | Shaft Cushion, Mobile 2 RPW |
| 28730 | V-Actor Set |
| 28740 | V-Actor V25 Transmitter |
| 28741 | V-Actor V35 Transmitter |
| 28745 | Carrying Bag |



